# Baby Care

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12

## Baby hygiene set with travel case

Please, read the instructions carefully and keep them for future reference.

### Features:

- Finger toothbrush.
- Nasal aspirator is designed to help relieve your baby's nasal congestion by brewing mucus from nose.
- The brush is made of soft nylon material.
- The comb has rounded tines, which provide an easy glide through the baby's hair.
- Scissors, Nail clipper, Nail file. Manicure Set was developed to offer comfort while grooming your child's nails.

### For your child's safety and health: WARNING!

#### • **Finger toothbrush**

**USE AND CARE:** Before first use sterilize at no more than 100° C. Put the finger toothbrush on your index finger. It must fit snugly onto the finger tip to make sure that it will not come out during use. Wet the bristles slightly with boiled warm water. Clean baby teeth, tongue and gums or massage gums using gentle pressure. Use for no more than 1 minute. After use wash in hot water with a mild detergent, rinse thoroughly, squeeze out and dry. Do not clean with solvents or harsh chemicals, do not wash in the dishwasher.

**WARNING!** This product must be handled only by an adult. The product is not a toy. Keep it out of reach of children. Check the product before each use. Throw away at the first signs of damage or weakness. Never use as a soother. Keep away from fire.

The product is intended for children 3 m+.

Content: Silicone

#### • **Nasal ASPIRATOR**

**USE AND CARE:** Before first use sterilize the silicone tip at no more than 100° C and wash the rubber pump with warm water and mild detergent. Rinse thoroughly, squeeze out and dry. Do not insert the silicone tip deep into the nasal cavity as this may cause injury. Squeeze the rubber pump and then gently insert the silicone tip in one of the baby's nostrils. Press the other nostril carefully with your fingertip. Gently reduce the pressure on the rubber pump, these will extract the nasal secretion from the baby's nose. After use wash in hot water with a mild detergent, rinse thoroughly, squeeze out and dry. Do not clean with solvents or harsh chemicals, do not wash in the dishwasher.

**WARNING!** This product must be handled only by an adult. The product is not a toy. Keep it out of reach of children. Check the product before each use. Throw away at the first signs of damage or weakness. Keep away from fire. For external use only.

The product is intended for children 0 m+.

Content: : PP, Silicone

#### • **Comb and Brush**

**USE AND CARE:** The comb and the hairbrush are intended for daily hair care of the baby. Keep the set clean, wash in warm water with a mild detergent. Rinse thoroughly and dry. Do not clean with solvents or harsh chemicals, do not wash in the dishwasher.

**WARNING!** This product must be handled only by an adult. The product is not a toy. Keep it out of reach of children. Check the product before each use. Throw away at the first signs of damage or weakness. Keep away from fire.

The product is intended for children 0 m+.

Content: PP,TPR

#### • **Scissors, Nail clipper, Nail file**

**USE AND CARE:** It is important to keep baby's nails clean and trimmed to prevent tearing and accidental scratching. Keep the nails straight and not too short. It is recommended to cut the nails while the baby is sleeping, when she/he is calm and does not move, to avoid the risk of injury. Keep the set clean, wash in warm water with a mild detergent. Rinse thoroughly and dry. Do not clean with solvents or harsh chemicals, do not wash in the dishwasher.

**WARNING!** This product must be handled only by an adult. The product is not a toy. Keep it out of reach of children. Check the product before each use. Throw away at the first signs of damage or weakness. Keep away from fire.

The product is intended for children 0 m+.

Content: PP,ABS, stainless steel

**WARNING!** Never leave the child unattended! These products should be used only by adults!

**Content:** PP, TPR, ABS, stainless steel, Silicone ; **Age:\*** Finger toothbrush: Age 3 m + / \* Nasal aspirator, hair and nail care sets: Age 0 m+ /

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